



Recombinant Ashbya gossypii RNA polymerase II transcription factor B subunit 5 (TFB5)

Product Code	CSB-YP762173DOT
Abbreviation	TFB5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q74ZQ0
Product Type	Recombinant Protein
Immunogen Species	Ashbya gossypii (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (Eremothecium gossypii)
Purity	≥85% (SDS-PAGE)
Sequence	MPRARKGALV QCDPSIRALI LQIDSGTHDI IWEELDETHL LVDPKVSYSI KEKLNWFLSK NIYNPSEEEEE NP
Source	Yeast
Target Names	TFB5
Protein Names	Recommended name: RNA polymerase II transcription factor B subunit 5 Alternative name(s): General transcription and DNA repair factor IIH subunit 5 Short name= TFIIH subunit 5
Expression Region	1-72
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.