



Recombinant *Ashbya gossypii* DNA polymerase epsilon subunit D (DPB4)

Product Code	CSB-BP762176DOT
Abbreviation	DPB4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q750A4
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MVCGESNAAA ISAKYVDNDP LQRPAARATP ANESHICDLM PPKGWKKDAQ GNYPTTSYIK EQEKLTVDDL LFPRSIITSL AKDAVHQAVQ TAEQDPRVML SKDASLALQR SSTVFNHLL MHARQIAQSN DRKSCSGEDV LKALDQIGLA GFESVVRERV VEYEKEVQRR RAEKTPAADE GQAEEGDAAD EEEGSHKRAK LDEH
Source	Baculovirus
Target Names	DPB4
Protein Names	Recommended name: DNA polymerase epsilon subunit D EC= 2.7.7.7 Alternative name(s): DNA polymerase II subunit D
Expression Region	1-204
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.