



# Recombinant *Ashbya gossypii* tRNA (adenine (58)-N (1))-methyltransferase catalytic subunit TRM61 (TRM61)

<b>Product Code</b>	CSB-EP762207DOT-B
<b>Abbreviation</b>	TRM61
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q755M8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) ( <i>Eremothecium gossypii</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSYFTEYKDL IAEGDLVLVW ISRGNVKQLR IKSGEVFNTR YGSFAHDDIV GRPYGSQIGV RTKGSNRFGE IHVLQPTPEL WTLSPHRTQ IVYTPDSSYI MQRLGCSPQS RVIEAGTGSG SFSHAFARTV GQLFSYEFHH MRYEQALAEF QEHGLIADGN VTITHRDVCQ KGFLIKAGDT TSYEFLEGQD SVSVNADCIF LDLPAPWEAI PHLDAVIAKD RTARLCCFSP CIEQVDKTL ALELHWEEI ETVEIQGRQY ESRRQMVRQL DDALERLRDV KRRKQHGTER RKRLGEQLGS VEELTDDVRP KTPKTSYNPF GKGSRVKEGD NGFEWKQVAK VESEIKSHTS YLTFASKILF KSRDENTVSQ LLEQYKDFNP NVKAASKLKD RPQETSFEES QASNSV
<b>Source</b>	<i>E. coli</i>
<b>Target Names</b>	TRM61
<b>Protein Names</b>	Recommended name: tRNA (adenine(58)-N(1))-methyltransferase catalytic subunit TRM61 EC= 2.1.1.220 Alternative name(s): tRNA(m1A58)-methyltransferase subunit TRM61 Short name= tRNA(m1A58)MTase subunit TRM61
<b>Expression Region</b>	1-406
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.