



# Recombinant *Ashbya gossypii* ATPase GET3 (GET3)

<b>Product Code</b>	CSB-MP762231DOT
<b>Abbreviation</b>	GET3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q759J2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) ( <i>Eremothecium gossypii</i> )
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MTDITPEASL RSLINSTTHK WIFVGGKGGV GKTTSSCSIA IQMALAQPTK QFLLISTDPA HNLSDAFNEK FGKDARKVTG MDNLSCMEID PSAALKDVND MAIANGGDDD GLSGLLQGGGA LADLTGSIPG IDEALSFMEV MKHIKKQEQG DGEHFDTVIF DTAPTGH TLR FLQLPTTLTK VLDKFGA IAG RLG PMLNSFA GNPNVDVLGK MNELKESVQK IKKQFTDPDL TTFVCVCISE FLSLYETERL IQELISYDMD VNSIIVNQLL FAESDKEHNC KRCQARWKMQ KKYLSQIDEL YEDFHIVKMP LCAGEIRGLE NLKKFSCFLN NKYDPETDGE LUYQLEQGL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	GET3
<b>Protein Names</b>	Recommended name: ATPase GET3 EC= 3.6.-.- Alternative name(s): Arsenical pump-driving ATPase Arsenite-stimulated ATPase Golgi to ER traffic protein 3 Guided entry of tail-anchored proteins 3
<b>Expression Region</b>	1-349
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.