



Recombinant *Ashbya gossypii* ATPase GET3 (GET3)

Product Code	CSB-YP762231DOT
Abbreviation	GET3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q759J2
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MTDITPEASL RSLINSTTHK WIFVGGKGGV GKTTSSCSIA IQMALAQPTK QFLLISTDPA HNLSDAFNEK FGKDARKVTG MDNLSCMEID PSAALKDVND MAIANGGDDD GLSGLLQGGA LADLTGSIPG IDEALSFMEV MKHIKKQEQG DGEHFDTVIF DTAPTGHTLR FLQLPTTLTK VLDKFGAIAG RLGPMNLNSFA GNPNVDVLGK MNELKESVQK IKKQFTDPDL TTFVCVCISE FLSLYETERL IQELISYDMD VNSIIVNQLL FAESDKEHNC KRCQARWKMQ KKYLSQIDEL YEDFHIVKMP LCAGEIRGLE NLKKFSCFLN NKYDPETDGE LUYQLEQGL
Source	Yeast
Target Names	GET3
Protein Names	Recommended name: ATPase GET3 EC= 3.6.-.- Alternative name(s): Arsenical pump-driving ATPase Arsenite-stimulated ATPase Golgi to ER traffic protein 3 Guided entry of tail-anchored proteins 3
Expression Region	1-349
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.