



Recombinant *Ashbya gossypii* Peptidyl-prolyl isomerase CWC27 (CWC27)

Product Code	CSB-YP762238DOT
Abbreviation	CWC27
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q75A74
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MSSEPSASGK CVLYTTKGEL QIELWAKECP KTVRSFLQSI QDGKWDGIVL GKVAEDAVWA PAGLCACSAE INGLRLRFNRR GLVGMAPGQD QPFFTLASRG ELDGRAVVFG TLVGQSVYRL MEIAQGEVGD DGKTFVYPAE VRRAEVTIPY FDGLSGQKRR AEPEQQAAPR PRKIATRVL EYEESEDEAS DPPLDVIRIRA AHDILQDERL TDALHGDSPV PRPAEAPRGS DSAPCVASPD RRAEAAPPTD GPPSRCASLA GPAPASPGPV TDLAPQLSAR EQGTLSSLAE FRCKRGGKNP LLG
Source	Yeast
Target Names	CWC27
Protein Names	Recommended name: Peptidyl-prolyl isomerase CWC27 Short name= PPIase CWC27 EC= 5.2.1.8 Alternative name(s): Rotamase CWC27
Expression Region	1-303
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.