



Recombinant Cysteine protease-like virA (virA)

Product Code	CSB-YP762461SZB
Abbreviation	virA
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q7BU69
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	≥85% (SDS-PAGE)
Sequence	<p>MQTSNITNHE RNDSSWMSTV KSTTEVSWNK LSFCDILLKI ITFGIYSPHE TLAEKHSEKK LMDSFSPSL S QDKMDGEFAH ANIDGISIRL CLNKGICSVF YLDGDKIQST QLSSKEYNNL LSSLPPKQFN LGKVHTITAP VSGNFKTHKP APEVIETAIN CCTSIIPNDD YFHVKDTDFN SVWHDYIRDI RASDSNSTKI YFNNIEIPLK LIADLINELG INEFIDSKKE LQMLSYNQVN KIINSNFPQQ DLCFQTEKLL FTSLFQDPAF ISALTSAFWQ SLHITSSSVE HIYAQIMSEN IENRLNFMPE QRVINNCGHI IKINAVVPMK DTAISASGGR AYEVSSSIPL SHITCNGVGI NKIETSYLVH AGTLPSSSEGL RNAIPPESRQ VSFAIISPDV</p>
Source	Yeast
Target Names	virA
Protein Names	Recommended name: Cysteine protease-like virA EC= 3.4.22.- Alternative name(s): Effector protein virA
Expression Region	1-400
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>