



Recombinant *Oryza sativa* subsp. japonica PHD finger protein ALFIN-LIKE 6 (Os01g0887700, LOC_Os01g66420)

Product Code	CSB-MP762512OFG
Abbreviation	Os01g0887700, LOC_Os01g66420
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7F2Z1
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MEGGGGGGGG GGGGGGGGGG GGAPYATRTA EEVFRDLRGR RAGMIKALTT DVEKFYKLCD PEKENLCLYG YPNETWEVTL PAEEVPPEIP EPALGINFAR DGMNEKDWLA LVAVHSDSWL LSVAFYFGAR FGFDREARRR LFNMNNLPT IFEVVTGAAK KQAKEKTPNS SSKSNKPSSK VQSKAESRSK SKLSAPKDEE GSGDDEGEEE EDDHDNTLCG TCGTNDGKDE FWICCDNCEK WYHGKCVKIT PARAEHIKQY KCPDCTNKRA RA
Source	Mammalian cell
Target Names	Os01g0887700
Protein Names	Recommended name: PHD finger protein ALFIN-LIKE 6
Expression Region	1-272
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.