



Recombinant *Drosophila melanogaster* Uncharacterized zinc finger protein CG12744 (CG12744)

Product Code	CSB-YP762551DLU
Abbreviation	CG12744
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7JUR5
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	METALAEVQL SCLLCEQTFD ATEKLDEHLP THFPQPVSTG QTCDICGRMT RSSLELHQHY KRYHEAHVPN TEGHFQCQLC DKVFLQLDHL KVHVKIEHAT DGYQPDEKSL DWQHYSRSL DTTNDYKLDP ILCPPPKRKY PPRSPFFNPN LWLGADNCFM
Source	Yeast
Target Names	CG12744
Protein Names	Recommended name: Uncharacterized zinc finger protein CG12744
Expression Region	1-160
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.