



Recombinant Drosophila melanogaster Ribosomal RNA-processing protein 8 (CG7137)

Product Code	CSB-MP762555DLU
Abbreviation	CG7137
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7K2B0
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MKPFEVPPWE EGADIEFEP VNPSGNGVAA KKKPKKKKPK KKKA AVPNEN LNVNKA AEKF TYRQGRGPLA PPQDSSDDDY EDVEDGIRAM HKVKSGRIEK QRPPTQATKG GK KELKLP E LAQA ALEAME AT SSTTPAAN SLASKLQSEL LGGRFRYINE QLYSTTSRKA EALFRKDSSA FEAYHAGYRQ QVEKWPINPL NRIIKTIKKI PKTAIIGDFG CGEGKLAQSV PNKVYSMDLV AARSDIACN ITDTPLQART LDVAVYCLSL MGTDLNEFFL EANRVLKLG TVYIAEIQSR FQDVREFVRC LNACGFDLNK KDVA VNYFYF FQFKKMRHVP KNTKMKAFSL KPCLYRKR
Source	Mammalian cell
Target Names	CG7137
Protein Names	Recommended name: Ribosomal RNA-processing protein 8 EC= 2.1.1.-
Expression Region	1-358
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.