



# Recombinant Human Serine/threonine-protein kinase MARK2 (MARK2)

<b>Product Code</b>	CSB-BP762568HU
<b>Abbreviation</b>	MARK2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7KZI7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSSARTPLPT LNERDTEQPT LGHLDSKPSS KSNMIRGRNS ATSADEQPHI GNYRLLKTIG KGNFAKVKLA RHILTGKEVA VKIIDKTQLN SSSLQKLFRE VRIMKVLNHP NIVKLFEVIE TEKTLYLVM E YASGGVEFDY LVAHGRMKEK EARAKFRQIV SAVQYCHQKF IVHRDLKAEN LLLDADMNIK IADFGFSNEF TFGNKLDTEFC GSPPYAAPEL FQGKKYDGPE VDVWSLGVIL YTLVSGSLPF DGQNLKELRE RVLRGKYRIP FYMSTDCENL LKKFLILNPS KRGTTLEQIMK DRWMNVGHED DELKPYVEPL PDYKDPRTTE LMVSMGYTRE EIQDSLVGQR YNEVMATYLL LGYKSSELEG DTITLKPRPS ADLTNSSAPS PSHKVQRSVS ANPKQRRFSD QAAGPAIPTS NSYSKKTQSN NAENKRPEED RESGRKASST AKVPASPLPG LERKKTTPPT STNSVLSTST NRSRNSPLLE RASLGQASIQ NGKDSLTPMG SRASTASASA AVSAARPRQH QKSMSASVHP NKASGLPTE SNCEVPRPST APQRVPVSP SAHNISSSGG APDRTNFPRG VSSRSTFHAG QLRQVRDQQN LPYGVTPASP SGHSQGRGA SGSIFSKFTS KFVRRNLSFR FARRNLNEPE SKDRVETLRP HVVGSGGNDK EKEEFREAKP RSLRFTWSMK TTSSMEPNEM MREIRKVLDA NSCQSELHEK YMLLCMHGTP GHEDFVQWEM EVCKLPRLSL NGVRFKRISG TSMFAKNIAS KIANELKL
<b>Source</b>	Baculovirus
<b>Target Names</b>	MARK2
<b>Protein Names</b>	Recommended name: Serine/threonine-protein kinase MARK2 EC= 2.7.11.1 EC= 2.7.11.26 Alternative name(s): ELKL motif kinase 1 Short name= EMK-1 MAP/microtubule affinity-regulating kinase 2 PAR1 homolog PAR1 homolog b
<b>Expression Region</b>	1-788
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a member of the Par-1 family of serine/threonine protein



kinases. The protein is an important regulator of cell polarity in epithelial and neuronal cells, and also controls the stability of microtubules through phosphorylation and inactivation of several microtubule-associated proteins. The protein localizes to cell membranes. Multiple transcript variants encoding different isoforms have been found for this gene.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.