



Recombinant Human Serine/threonine-protein kinase MARK2 (MARK2)

Product Code	CSB-EP762568HU-B
Abbreviation	MARK2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7KZI7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSSARTPLPT LNERDTEQPT LGHLDSKPSS KSNMIRGRNS ATSADEQPHI GNYRLLKTIG KGNFAKVKLA RHILTGKEVA VKIIDKTQLN SSSLQKLFRE VRIMKVLNHP NIVKLFEVIE TEKTLYLVM E YASGGVEFDY LVAHGRMKEK EARAKFRQIV SAVQYCHQKF IVHRDLKAEN LLLDADMNIK IADFGFSNEF TFGNKLDTF C GSPPYAAPEL FQGKKYDGPE VDVWSLGVIL YTLVSGSLPF DGQNLKELRE RVLRGKYRIP FYMSTDCENL LKKFLILNPS KRGTTLEQIMK DRWMNVGHED DELKPYVEPL PDYKDP RRE LMVSMGYTRE EIQDSL VGQR YNEVMATYLL LGYKSSELEG DTITLKPRPS ADLTNSSAPS PSHKVQRSVS ANPKQRRFSD QAAGPAIPTS NSYSKKTQSN NAENKRPEED RESGRKASST AKVPASPLPG LERKKTTP TP STNSVLSTST NRSRNSPLLE RASLGQASIQ NGKDSL TMPG SRASTASASA AVSAARPRQH QKSMSASVHP NKASGLP PTE SNCEVPRPST APQRVPV ASP SAHNISSSGG APDR TNFPRG VSSRSTFHAG QLRQVRDQQN LPYGVTPASP SGHSQGRRGA SGSIFSKFTS KFVRRNLSFR FARRNLNEPE SKDRVETLRP HVVGSGGNDK EKEEFREAKP RSLRFTWSMK TTSSMEPNEM MREIRKVLDA NSCQSELHEK YMLLCMHGTP GHEDFVQWEM EVCKLPRLSL NGVRFKRISG TSM AFKNIAS KIANELKL
Source	E.coli
Target Names	MARK2
Protein Names	Recommended name: Serine/threonine-protein kinase MARK2 EC= 2.7.11.1 EC= 2.7.11.26 Alternative name(s): ELKL motif kinase 1 Short name= EMK-1 MAP/microtubule affinity-regulating kinase 2 PAR1 homolog PAR1 homolog b
Expression Region	1-788
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a member of the Par-1 family of serine/threonine protein



kinases. The protein is an important regulator of cell polarity in epithelial and neuronal cells, and also controls the stability of microtubules through phosphorylation and inactivation of several microtubule-associated proteins. The protein localizes to cell membranes. Multiple transcript variants encoding different isoforms have been found for this gene.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.