



# Recombinant Human MOB kinase activator 1B (MOB1B)

<b>Product Code</b>	CSB-EP762574HU
<b>Abbreviation</b>	MOB1B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7L9L4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SFLFGSRSS KTFKPKKNIP EGSHQYELLK HAEATLGSGN LRMAVMLPEG EDLNEWVAVN TVDFFNQINM LYGTITDFCT EESCPVMSAG PKYEYHWADG TNIKKPIKCS APKYIDYLMT WVQDQLDDET LFPSKIGVPF PKNFMSVAKT ILKRLFRVYA HIYHQHFDPV IQLQEEAHLN TSFKHFIFFV QEFNLIDRRE LAPLQELIEK LTSKDR
<b>Source</b>	E.coli
<b>Target Names</b>	MOB1B
<b>Protein Names</b>	Recommended name: MOB kinase activator 1B Alternative name(s): Mob1 homolog 1A Short name= Mob1A Mob1B Mps one binder kinase activator-like 1A
<b>Expression Region</b>	2-216
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein is similar to the yeast Mob1 protein. Yeast Mob1 binds Mps1p, a protein kinase essential for spindle pole body duplication and mitotic checkpoint regulation.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.