



Recombinant Porphyromonas gingivalis Ribosome-recycling factor (frr)

Product Code	CSB-MP762745EYA
Abbreviation	frr
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7MTP6
Product Type	Recombinant Protein
Immunogen Species	Porphyromonas gingivalis (strain ATCC BAA-308 / W83)
Purity	≥85% (SDS-PAGE)
Sequence	MEIKELIAKA TANMEKAVEY LDEQLSHVRA GKASPKLLDG IMVMYYGNAT PLTQVASINT PDAKTIVVTP WERSLIKDIE KAIMDSPLGI TPENNGELIR LGLPPLTEER RRQLVKQTKG DAEDAKVSVR NARRDAIDAI KKS VKTDGTP EDVAKDAEAE MQKVH DRYIK KIDELFAEKE KEIMTV
Source	Mammalian cell
Target Names	frr
Protein Names	Recommended name: Ribosome-recycling factor Short name= RRF Alternative name(s): Ribosome-releasing factor
Expression Region	1-186
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.