



Recombinant Photorhabdus luminescens subsp. laumondii DnaA-homolog protein hda (hda)

Product Code	CSB-YP762802PIJ
Abbreviation	hda
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7N3G5
Product Type	Recombinant Protein
Immunogen Species	Photorhabdus luminescens subsp. laumondii (strain DSM 15139 / CIP 105565 / TT01)
Purity	≥85% (SDS-PAGE)
Sequence	MNTPSQLSLP LYLPDDETF SFFPGENATL LAAIKLAINQ PHGSYIYFWS RDSGGRSHLL HAACAELSLK DEAVGYVPLD KRAYFIPEVL DGMEHLSLVC IDNIESIAGD EEWEMAIFNL YNRILEIGRT CLLISGDRPP RQINLKLPDL ASRLDWGQIY KLQPLSDDEK LQALQLRAKL RGFELPEDVG RFLKRLDRE MRTLFMALNQ LDHASIVAQR KLTIPFVKDI LHL
Source	Yeast
Target Names	hda
Protein Names	Recommended name: DnaA-homolog protein hda
Expression Region	1-233
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.