



# Recombinant Photorhabdus luminescens subsp. laumondii Siroheme synthase (cysG)

<b>Product Code</b>	CSB-YP762831PIJ
<b>Abbreviation</b>	cysG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7N8L2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Photorhabdus luminescens subsp. laumondii (strain DSM 15139 / CIP 105565 / TT01)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDYLPFVVEL KGRLVLLVGG GEVAARKATL LLRAGALLQV VAPELCSELQ QRYQAGELEW YQGEFQPEYL DGIFLVIAAT DDRILNHQVF SEADKRSILV NVVDDQVHCS FIFPSIDRS PVLVAISSAG KAPVLARLIR EKLEALLPSS LGTMAKIAGK WRERVKQRLT SMRQRSSFWE QAFNGRFAML VANGQIQQAE KQLEQQLEQS DLQGELALVG AGPGDPGLLT LKGLQVIQQA DVVLYDHLVS SDVLDLIRRD ADKICVGKRA GNHSVSQEET NRLIVKFARQ GKKVVRLKGG DPFIFGRGGE ELQVAASAGI PFQVVPGITA AIGATAYAGI PLTHREHSQS ITFITGHCRE KGNELDWPAL ARGHQTLVIY MGTVKAALIS HQLILHGRAE DTPVAVIGCG TRLEQQVLTG TLLELEQLAQ QAPSPALLVV GEVAQLHHQI AWFGQQSIAK ISRPAVVDF A
<b>Source</b>	Yeast
<b>Target Names</b>	cysG
<b>Protein Names</b>	Recommended name: Siroheme synthase Including the following 3 domains: Uroporphyrinogen-III C-methyltransferase Short name= Urogen III methylase EC= 2.1.1.107 Alternative name(s): SUMT Uroporphyrinogen III methylase S
<b>Expression Region</b>	1-470
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.