



Recombinant Mycoplasma gallisepticum Peptide deformylase (def)

Product Code	CSB-EP762849MWI
Abbreviation	def
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7NAK8
Product Type	Recombinant Protein
Immunogen Species	Mycoplasma gallisepticum (strain R(low / passage 15 / clone 2))
Purity	>85% (SDS-PAGE)
Sequence	MKSKLKPTND WLVTDDNPKM REVCTEVKFP LSQEVLDIID KMLAYVDESFD DDNAEKYDIR PGIGIAANQL GLNQRFFYVH FTDFCQKEHR YLLINPEWID KSLNKAYLAV GEGCLSVPKD KDGIVIRSET VKLKGFDYLT QKDVEISAHG LLAMCLQHEM DHLEGKFYYD SINMMKPYFK KDEWVSIDQK VCDQCK
Source	E.coli
Target Names	def
Protein Names	Recommended name: Peptide deformylase Short name= PDF EC= 3.5.1.88 Alternative name(s): Polypeptide deformylase
Expression Region	1-196
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.