



Recombinant *Chromobacterium violaceum* Transcription termination factor Rho (rho)

Product Code	CSB-MP762957CKA
Abbreviation	rho
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7NXP1
Product Type	Recombinant Protein
Immunogen Species	<i>Chromobacterium violaceum</i> (strain ATCC 12472 / DSM 30191 / JCM 1249 / NBRC 12614 / NCIMB 9131 / NCTC 9757)
Purity	≥85% (SDS-PAGE)
Sequence	MHLSDLKHLH VSELVEMAIS NEIDGANRLR KQDLIFALLK NQAKKGESIF GEGTLEVLPD GFGFLRSPDT SYLAGPDDIY VSPSQIRRFN LHTGDSIEGE IRTPKDGERY FALVKVDKVN GEAPENSKHK ILFENLTPLF PTEQFKLERE IRAEENITGR IIDLISPIGK QQRALLVAPP KSGKTVMLQH IAHAITANHP EAVLIVLLID ERPEEVTEMQ RSVRGEVVSS TFDEPATRHV QVAEMVIEKA KRLVEHKKDV VILLDSITRL ARAYNTVVPA SGKVLTTGGVD ANALQRPKRF FGAARNVEEG GSLTIVATAL IDTGSRMDDV IYEEFKGTGN SEIHLDRRMA EKRIFFPALNI NRSCTRREEL LVPQDQLQRI WVLRKLLYPM DDLEAMEFLQ DKIKATKSNQ AFFDSMRR
Source	Mammalian cell
Target Names	rho
Protein Names	Recommended name: Transcription termination factor Rho EC= 3.6.4.- Alternative name(s): ATP-dependent helicase Rho
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.