



Recombinant Chromobacterium violaceum Ribosomal RNA small subunit methyltransferase A (rsmA)

Product Code	CSB-YP762975CKA
Abbreviation	rsmA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7P1U1
Product Type	Recombinant Protein
Immunogen Species	Chromobacterium violaceum (strain ATCC 12472 / DSM 30191 / JCM 1249 / NBRC 12614 / NCIMB 9131 / NCTC 9757)
Purity	>85% (SDS-PAGE)
Sequence	MSKHIPRKRF GQNFLQDASV IAGIVHAVNP QPDDIVIEIG PGLGAIKPL LARLKHHLHV EIDRDIERL KAEHPADKLT IHAGDALAFD FASVSEAPLK IVGNLPYNIS TPLLHFLASY GNRVTDMHFM LQKEVIERMV AEPSTADYGR LTVMLQYRFY MENILFVPE AFWPPPKVDS AVVRMIPAPG RCGTARDEAL LEKLVSQAFA QRRKTLRNNL KGLADAADLE ALGIDPGLRP ENLPVEDFVR LANHLHDKGA RG
Source	Yeast
Target Names	rsmA
Protein Names	Recommended name: Ribosomal RNA small subunit methyltransferase A EC= 2.1.1.182 Alternative name(s): 16S rRNA (adenine(1518)-N(6)/adenine(1519)-N(6))-dimethyltransferase 16S rRNA dimethyladenosine transferase 16S rRNA dimethylase
Expression Region	1-262
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.