



Recombinant Danio rerio UPF0489 protein C5orf22 homolog (zgc:56556)

Product Code	CSB-EP763067DIL
Abbreviation	zgc:56556
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7SZF1
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MNAPPQKR VY PKLPVWIVED HHEVVQHIYR AIGSRHIPMK GIKMVHLD SH PDLLIPVNMP ADTVYDKETL LRELSIENWI MPMVYAGHVS HVAWLHPYWA QQIREGQHCM CVGKDSSTTT IRVTSKDDYF LSDALYVPLD HLEHPKELHL HVIRVDPVVS SQTTKLENGQ SGAKIPKAAQ TQDDMQSKAD TPCTSSSQPP DGSAASGNIS ETAKKKADDG STSYITDKLL AVIEQTDPFI LDIDLDFESC KNPFKDMFSQ EEFTLLKELY SFSKPQQDPD EEELLECVER RTRQLEDLEA AFADLLEDDG QETVERLAAN PGMKSLFRLV HSLKNRPSPP DYEMVHQAGL TCDNSELP HH ISSDEEIQQM ITAVQLFLET LPKPTIVTIS RSSLDEYCPA EQVDSIQIRV LDILESLFGC LDVHRDYESI PAESTSHTA
Source	E.coli
Target Names	zgc:56556
Protein Names	Recommended name: UPF0489 protein C5orf22 homolog
Expression Region	1-439
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.