



# Recombinant Danio rerio ERO1-like protein alpha (ero1l)

<b>Product Code</b>	CSB-YP763085DIL
<b>Abbreviation</b>	ero1l
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7T3D1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	GSAAHRCFCQ VTGTLDDCAC DVETIDKFNN KDIFPKLQKL LSSDYFRFYK VNLNNGCPFW TDHSQCGLKY CAVKPCSPDE VPEGLKSSSY KYSEKASHDT EECEKAEKLG AVNGSLSDET RQALEEWKKY DDESDRFCML DDEDSPESQY VDLLLLNPERF TGYKGAEAWR IWNSIYEENC FKPYSVNRPL NPLASNSGDD GQGFYRWLEG LCVEKRAFFR LISGLHASIN IHLRALLD ENWFEMKWGH NVSEFQQRFD EDLTKEGEPK RLRNLYFLYL IELRALAKIL PYFERSTFQL YTGQDTQDDQ NKKLLELLH VAKSFPLHFD ETALFAGNNK EAMKLKEDFK LTFKNISRIM DCVECFKCRL WGKLQTQGLG TALKILFSEK QIEAMPNKQH QSIISAQSAG NCVFVQRFRK NLHKCQRVEE LQIVVETQT VTKMFMQERA TVQTDMMHKHM TFPFMFGYIQ
<b>Source</b>	Yeast
<b>Target Names</b>	ero1a
<b>Protein Names</b>	Recommended name: ERO1-like protein alpha Short name= ERO1-L Short name= ERO1-L-alpha EC= 1.8.4.- Alternative name(s): Endoplasmic oxidoreductin-1-like protein Oxidoreductin-1-L-alpha
<b>Expression Region</b>	21-489
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.