



# Recombinant Prochlorococcus marinus subsp. pastoris GTPase Der (der)

<b>Product Code</b>	CSB-EP763367EYQ
<b>Abbreviation</b>	der
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7V2S6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Prochlorococcus marinus subsp. pastoris (strain CCMP1986 / NIES-2087 / MED4)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MTLPSIAIIG RPNVGKSTLV NRLCQSNDAI VFDKPGVTRD RTYQNASWAG REFQVVDTGG LVFEDDSEFL PEIRTQVFLA LEEASLALFV VDGNGQVTDG DLSIAKWLRN SACKTIVAVN KCESTSLGVS LASEFWKLGL GEPYPVSAIH GSGTGDLLDL VIDGFPKDLN VEDKEDKIMM SIIGRPNVGK SLLNAICGE KRAIVSDISG TTTDSIDTLI KKDSHLWKIV DTAGIRRKKN VKYGTEFFGI NRAFKSIDRS DVCVLVIDAI DGVTDDQDQKL AGRIEEQGRA CVIVVNKWDL VEKNNSTIYQ VEKELRSKLY FLHWSKMIFI SALTGQRVQN IFEHALSAVT QHRRRVTTTSV VNEVLKEALG WKSPPTKRSG KQGRLYYGTQ VKNQPPTFTL FVNDPKLFGI TYRRYIEKQI RLNLGFEGSP IILLWRGKQQ RDLEKETSCK NINIIQKD
<b>Source</b>	E.coli
<b>Target Names</b>	der
<b>Protein Names</b>	Recommended name: GTPase Der Alternative name(s): GTP-binding protein EngA
<b>Expression Region</b>	1-458
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.