



Recombinant Candida glabrata Enolase 1 (ENO1)

Product Code	CSB-MP763769CZI
Abbreviation	ENO1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6FTW6
Product Type	Recombinant Protein
Immunogen Species	Candida glabrata (strain ATCC 2001 / CBS 138 / JCM 3761 / NBRC 0622 / NRRL Y-65) (Yeast) (Torulopsis glabrata)
Purity	>85% (SDS-PAGE)
Sequence	MSVTKVHART VYDSRGNPTV EVEVYTKDGM FRAIVPSGAS TGSHEALELR DGDKSKWEGK GVTKAVSNVN DIIGPKLIES KLDVKDQKAI DEFMIKLDGT KNKSKLGANA ILGVSLAVAR AGAAAKGVPL YQHIAELADM KKEPYVIPCP FFNVLNGGVH AGGNLALQEF LIAPVGAESF HEALRLGSEV YHKLKALAKK RIGSSAGNVG DEGGIAPSL S TPFEALDIY DAIKEAGHEG KVKIAMPAS SEFFQGDKYD LDFKNPHPDA KNKLSGAQLG DYYKTILEKY PIVSLEDPFA EDDWEAWTNF FPKAGVQIIA DDLTVTNPER IQTAIDKKA DCLLLKVNQI GSLTESINSA KLAYGAGWGV QVSHRSGETE DTFIADLVVG LRTGQIKSGS LARSERLAKW NQILRIEEEL GASNTIFAGA KFHKGQLL
Source	Mammalian cell
Target Names	ENO1
Protein Names	Recommended name: Enolase 1 EC= 4.2.1.11 Alternative name(s): 2-phospho-D-glycerate hydro-lyase 1 2-phosphoglycerate dehydratase 1
Expression Region	1-438
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.