



Recombinant *Xenopus tropicalis* ELAV-like protein 1 (elavl1)

Product Code	CSB-MP764089XBF
Abbreviation	elavl1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6GLB5
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus tropicalis</i> (Western clawed frog) (<i>Silurana tropicalis</i>)
Purity	>85% (SDS-PAGE)
Sequence	MSNGYGDHMD DVCRRDDIGRT NLIVNYLPQN MTQDELRS LF SSIGEVESAK LIRDKVAGHS LYG FVNYLN AKDAERAIN T LNGLRLQSKT IKVSVAR PSS ESIKDANLYI SGLPRTMTQK DVEDMFLPFG RIINSRVLVD QATGLSRGVA FIRFDKRSEA EEAIASFNGH KPPGSSEPIT VKFAANPNQN KNMALLSQLC HSPARRFGGP VHHQAQRFRF SPMGVDHMSS ISGVNVASSA SSGWCIFIYN LGQDADEGIL WQMFGPFGAV TNVKVIRDFN TNKCKGFGFV TMTNYEEAAM AIASLNGYRL GDKTLQVFFK TSKSHK
Source	Mammalian cell
Target Names	elavl1
Protein Names	Recommended name: ELAV-like protein 1 Alternative name(s): Protein ElrA
Expression Region	1-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.