



# Recombinant Rat Keratin, type I cytoskeletal 24 (Krt24)

<b>Product Code</b>	CSB-MP764235RA
<b>Abbreviation</b>	Krt24
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6IFX1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MFCSAQKGGSC SSRVSSSGAV GSRGCTGSGS SYGLGGGSAW GFQSSSSSWG LSGGSKGSIG GGFSSCSVRG GFGAGSSFGG GSGFGGGSSG GVSSYGGSLG GGLGGIGGYD GLLSGSEKQ TMQGLNDRLA NYLDKVRAL EANTDLET KI KDWYGRHGSG KDGPGRDYSQ YCSVIEDLKN QIISATCENA RLALQIDNAR LAADDFRMKY EHELCLRQSL EADINGLRKV LDEMTMTRCD LEMQIEGLTE ELVFLRKNHE EEMKCLQGSS GGDVTVEMNA TPGTDLTKLL NDMRAQYEAM AEQNRQEA EK QFNERSACLQ AQISTDAGAA NCAKSEVMEL RRTVQTLEIE LQSQLALKCS LEGLADTEA RYVAQLSGIQ TQISSLEEQL SQIRGETQCQ SAEYECLLDI KTRLEQEIET YRLLNGDGG GCDYRNLVSR QVVLNDSNFG SCSGQEKDPS KTRVTKTIIE EVVDGKVVSS QVSNISEVKI K
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Krt24
<b>Protein Names</b>	Recommended name: Keratin, type I cytoskeletal 24 Alternative name(s): Cytokeratin-24 Short name= CK-24 Keratin-24 Short name= K24 Type I keratin-24
<b>Expression Region</b>	1-491
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a member of the type I (acidic) keratin family, which belongs to the superfamily of intermediate filament (IF) proteins. Keratins are heteropolymeric structural proteins which form the intermediate filament. These filaments, along with actin microfilaments and microtubules, compose the cytoskeleton of epithelial cells. The type I keratin genes are clustered in a region of chromosome 17q12-q21.



---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.