



Recombinant Danio rerio Phosphatidylinositol 5-phosphate 4-kinase type-2 gamma (pip4k2c)

Product Code	CSB-YP764249DIL
Abbreviation	pip4k2c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6IQE1
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MASLGNSGSA SSPMVMLAPK TKTKKRHFMQ QKVKVFASD PMLSVFMWGV NHSINDLNQV PVPVMLLPDD FKANTKIKVN NHLFNKENLP SHFEFKEYCP QVFRNLRERF GIEDLDYQAS LARSAPMKGD GQGEGLLFTS YDRTLIVKQI SSEEVADMHN ILSEYHQHIV KCHGSTLLPQ FLGMYRITVE SEDTYLIVMR NMFSHRLLVH RKYDLKGSV DREASDKEKV KELPTFKDMD FRNNMQKVYV TEEQKEKMM E KLN RDVEFLV KLKIMDYSLL LGIHDVARGE REEEAEPC YEDDADPENG LAPALQVGSY GTSPEGIAGY MNSIKPLGPG EFDPYIDVYA VKSAPGAPQR EVYFMGLIDV LTQYDTKKKA AHAAKTVKHG AGAEISTVHP EQYAKRFREF ISNIFA
Source	Yeast
Target Names	pip4k2c
Protein Names	Recommended name: Phosphatidylinositol 5-phosphate 4-kinase type-2 gamma EC= 2.7.1.149 Alternative name(s): Phosphatidylinositol 5-phosphate 4-kinase type II gamma Short name= PI(5)P 4-kinase type II gamma Short name= PIP4KII-
Expression Region	1-416
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.