



# Recombinant Mouse Multivesicular body subunit 12B (Fam125b)

<b>Product Code</b>	CSB-EP764286MO-B
<b>Abbreviation</b>	Fam125b
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6KAU4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MRSCFCVRRS RDP PPPQPPP PQRGTDQATM PEVKELSEAL PETPMDPITG VGVVASRNRA PTGYDVVAQT ADGVDADLWK DGLFKSKVTR YLCFTRSF ENSHLGNLV DMKLIDVKDT LPVGFPIQE TVDTQEVVFR KKRLCIKFIP RDSTEEAICD IRIMGR TKQA PPQYTFIGEL NSMGIWYRMG RVPRNHDS SQ PTTSPQSSAS STPAPNLPRH ISLTLPATFR GRNNTSTDYE YQLSNLY AIS AMDGVPFMIS EKFSICIPESM QPFDLLGITI KSLAEIEKEY EYSFRTE QSA AARLPPSPTR CQQIPQS
<b>Source</b>	E.coli
<b>Target Names</b>	Mvb12b
<b>Protein Names</b>	Recommended name: Multivesicular body subunit 12B Alternative name(s): ESCRT-I complex subunit MVB12B Protein FAM125B
<b>Expression Region</b>	1-317
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.