



Recombinant Mouse Multivesicular body subunit 12B (Fam125b)

Product Code	CSB-MP764286MO
Abbreviation	Fam125b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6KAU4
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MRSCFCVRRS RDP PPPQPPP PQRGTDQATM PEVKELSEAL PETPMDPITG VGVVASRNRA PTGYDVVAQT ADGVDADLWK DGLFKSKVTR YLCFTRSFSK ENSHLGNVLV DMKLIDVKDT LPVGFPIQEQ TVDTQEVVFR KKRLCIKFIP RDSTEEAICD IRIMGR TKQA PPQYTFIGEL NSMGIWYRMG RVPRNHDSSQ PTTPSQSSAS STPAPNLPRH ISLTLPATFR GRNNTSTDYE YQLSNLYAIS AMDGVPFMIS EKFCIPESM QPFDLLGITI KSLAEIEKEY EYSFRTEQSA AARLPPSPTR CQQIPQS
Source	Mammalian cell
Target Names	Mvb12b
Protein Names	Recommended name: Multivesicular body subunit 12B Alternative name(s): ESCRT-I complex subunit MVB12B Protein FAM125B
Expression Region	1-317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.