



Recombinant Arabidopsis thaliana F-box/FBD/LRR-repeat protein At2g04230 (At2g04230)

Product Code	CSB-YP764642DOA
Abbreviation	At2g04230
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NKX3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MEQKFKTDSM NEDRISDLPD ALLLQILSSL PTENAIATSV LSKRWRSWLT MLPKLKFDSN FNPVFDDNI DPTMFSENVY KTLSLHKAPV LESLHLSFEG RTDCLHVGIW IATAFARGVR KVLVDSFYQE DQVTPLPSVL FSYNDSLEIL KPKCAIDLDF PSRVCLKSLR KLYLDQVHFK DEESVCNLLC GCPSLQDLVV HRYSNADVAT FTIASPSLQR LTIEDLRQEG GYNGSYVIN APGLKYLIN GVIDIESCLI DKALELEEAK ISNVSGITNE NILESLSAK RLILHLSPLE VKVPTGKIFD QLGCELLTH EREWVLLSI MLDSSPKLQI LKLTDVYLHD NKTNPDERKW NPPKCAPECL LFHLETFLWI GYEWQRGDEK EVATYILENA RRLKKATFST KRIGREKLED FEKRREMLNE LAIVLWDSNS CHLVFEST
Source	Yeast
Target Names	At2g04230
Protein Names	Recommended name: F-box/FBD/LRR-repeat protein At2g04230
Expression Region	1-448
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.