



Recombinant Arabidopsis thaliana Protein ABIL3 (ABIL3)

Product Code	CSB-EP764646DOA
Abbreviation	ABIL3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NMC6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSAAATMPMP REASNYDEIS MQQSMLFSDS LKDLKNLRTQ LYSAAEYFEL SYTNDEQKQI VVETLKDYAI KALVNTVDHL GSVTYKVNDV VDEKVDEVAG TELRVSCIEQ RLRMCQEYMD HEGRSQQSLV IDTPKFHKRY FLPSGEIKRG GNLAKLKNVE GSFDGEDDWN QFRNAVRATI RETPPPPVRK PILQSPSQRK PQRSATFSFS SIATAPKKEQ DKRAVSPHRF PLLRSGSVAI RPSSISRPTT PSKSRVAVTPT PKRYPSEPRR SASVRVAFEK EAQKEPEHQQ QPSKSKRLLK ALLSRRKTKK DDTLYTYLDE Y
Source	E.coli
Target Names	ABIL3
Protein Names	Recommended name: Protein ABIL3 Alternative name(s): Abl interactor-like protein 3 Short name= AtABIL3
Expression Region	1-321
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.