



Recombinant Arabidopsis thaliana F-box/kelch-repeat protein At5g26960 (At5g26960)

Product Code	CSB-EP764649DOA
Abbreviation	At5g26960
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NPN5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSENCNSRHF SWLMKSCLPN PSDAKSLVQI HQPSSTAANS SATIASLPDD LLEICISRVP SSSIPSLAVV CRRWSRLLHS PYFLHLRRRL GLLRHSLFAI STVDSGLFAA DLQFQSEIAS WKVSLAVSSR SVGVDGSYGS LSHARAAAIG PRVYVVS RNA V LRYDSWMGT LNLRSPIFIP RKKFAIAVVS GKIYVAGGGG GSEVAAAVEE YDPELNRWEV VTQSARKRYG CIGAAVDGVF YVIGGLKIGN ETRAVAARA YASSMDLFDV ESRQWLSRS VPGGGCVVAA CAAVGYVYVL TSHAVELSFV RFDARRRGGN SGFGEWQRLK SPPLPAQVRL DGTVRFSCVG VEDKVAVVQV VGCIDLLRR SGRGERGIRE SLVLLYDTTE GEWRRAADLP EMITRAACAC VEW
Source	E.coli
Target Names	At5g26960
Protein Names	Recommended name: F-box/kelch-repeat protein At5g26960
Expression Region	1-413
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.