



# Recombinant Arabidopsis thaliana F-box/kelch-repeat protein At5g26960 (At5g26960)

<b>Product Code</b>	CSB-YP764649DOA
<b>Abbreviation</b>	At5g26960
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6NPN5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSENCNSRHF SWLMKSCLPN PSDAKSLVQI HQPSSTAANS SATIASLPDD LLEICISRVP SSSIPSLAVV CRRWSRLLHS PYFLHLRRRL GLLRHSLFAI STVDSGLFAA DLQFQSEIAS WKVSLAVSSR SVGVDGSYGS LSHARAAAIG PRVYVVS RNA VLR YDSW MGT LNL RSP MIFP RKKFAIAVVS GKIYVAGGGG GSEVAAAVEE YDPELNRWEV VTQSARKRYG CIGAAVDGVF YVIGGLKIGN ETRAVAARA YASSMDLFDV ESRQWLSRS VPGGGCVVAA CAAVGYVYVL TSHAVELSFWRFDARRRGGN SGFGEWQRLK SPPLPAQVRL DGTVRFSCVG VEDKVAVVQV VGCIDLLRR SGRGERGIRE SLVLLYDTTE GEWRRAADLP EMITRAACAC VEW
<b>Source</b>	Yeast
<b>Target Names</b>	At5g26960
<b>Protein Names</b>	Recommended name: F-box/kelch-repeat protein At5g26960
<b>Expression Region</b>	1-413
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.