



Recombinant Danio rerio Methionine-R-sulfoxide reductase B2, mitochondrial (msrb2)

Product Code	CSB-MP764677DIL
Abbreviation	msrb2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NW52
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	SLTRYDQST DWQRKLSPEQ YVVTREKGT VPFSGIYLNH NEVGMVYHCVC CDSPLFSSEA KYDAGAGWPS FCEAHGTWEK DESHANIVRR PDNSLGSTGT EVICKHCDAAH LGHVFEDGPD PTGQRFCINS VALNFKPREK
Source	Mammalian cell
Target Names	msrb2
Protein Names	Recommended name: Methionine-R-sulfoxide reductase B2, mitochondrial Short name= MsrB2 EC= 1.8.4.-
Expression Region	42-180
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.