



Recombinant Human TELO2-interacting protein 2 (TTI2)

Product Code	CSB-EP764687HU
Abbreviation	TTI2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NXR4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MELDSALEAP SQEDSNLSEE LSHSAFGQAF SKILHCLARP EARRGNVKDA VLKDLGLDIE ATEFDRLFEG TGARLRGMPE TLGQVAKALE KYAAPSKEEE GGGDGHSEAA EKAAQVGLLF LKLLGKVETA KNSLVGPAWQ TGLHHLAGPV YIFAIHSLE QPWTTPRSRE VAREVLTSLL QVTECGSVAG FLHGENEDEK GRLSVILGLL KPDLYKESWK NNPAIKHVFS WTLQQVTRPW LSQHLERVLP ASLVISDDYQ TENKILGVHC LHHIVLNVPA ADLLQYNRAQ VLYHAISNHL YTPEHHLIQA VLLCLLDLFP ILEKTLHWKG DGARPTTHCD EVLRILITHM EPEHRLLLRR TYARNLPAPV NRLGILTVRH LKRLERVIIG YLEVYDGPEE EARLKILETL KLLMQHTWPR VSCRLVLLK ALLKLICDVA RDPNLTPESV KSALLQEATD CLILLDRCSQ GRVKGLLAKI PQSCEDRKVV NYIRKVQQVS EGAPYNGT
Source	E.coli
Target Names	TTI2
Protein Names	Recommended name: TELO2-interacting protein 2
Expression Region	1-508
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.