



# Recombinant Human TELO2-interacting protein 2 (TTI2)

<b>Product Code</b>	CSB-MP764687HU
<b>Abbreviation</b>	TTI2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6NXR4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MELDSALEAP SQEDSNLSEE LSHSAFGQAF SKILHCLARP EARRGNVKDA VLKDLGLDIE ATEFDRLFEG TGARLRGMPE TLGQVAKALE KYAAPSKEEE GGGDGHSEAA EKAAQVGLLF LKLLGKVETA KNSLVGPAWQ TGLHHLAGPV YIFAIHSLE QPWTTPRSRE VAREVLTSLL QVTECGSVAG FLHGENEDEK GRLSVILGLL KPDLYKESWK NNPAIKHVFS WTLQQVTRPW LSQHLERVLP ASLVISDDYQ TENKILGVHC LHHIVLNVPA ADLLQYNRAQ VLYHAISNHL YTPEHHLIQA VLLCLLDLFP ILEKTLHWKG DGARPTTHCD EVLRILITHM EPEHRLRR TYARNLPAFV NRLGILTVRH LKRLERVIIG YLEVYDGPEE EARLKILETL KLLMQHTWPR VSCRLVLLK ALLKICDVA RDPNLTPESV KSALLQEATD CLILLDRCSQ GRVKGLLAKI PQSCEDRKVV NYIRKVQQVS EGAPYNGT
<b>Source</b>	Mammalian cell
<b>Target Names</b>	TTI2
<b>Protein Names</b>	Recommended name: TELO2-interacting protein 2
<b>Expression Region</b>	1-508
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.