



Recombinant Danio rerio Interferon regulatory factor 2-binding protein 2-A (irf2bp2a)

Product Code	CSB-EP764702DIL-B
Abbreviation	irf2bp2a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NZT6
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSSAAVAASR RQSCYLCDLP RMPWAMIWDF TEPVCRGCVN YEGADRIEFV IETARQLKRA HGFQDGRPAG PGKQLPGKEI HASGDPGSRP PQPLDRYPLS SERPPRLGAE YQSARQANGI PVPNGFSKPD DPPELNRQSP NPRRTSGVPP NLVPLVNGSM AHAMNGRPGQ MSDLGGNKRP ASVSSTEHDK DKHRPDSLTP EISDGHKTRA DDWLNKNKTV RELMTLNPHYD GRFKKEHMQQ RVMYETGPAS KSDRGKHPKN LKRKASPEPD GEGTAGKING DGQQWLPSAS GDGLKMPPVP PPSFTSPPST ISPHPRTPP EAGAAAAATQ NGHSPIAALI LAADNAGGNS SPK DANQVHS TTRRNSNSPL SPNQRRLAQG AGTPHVPGMD AVHPQSIPDS SVPSSIPLCC TLCHERLEDT HFVQCPSVPS HKFCFPCSRE SIKQQGATGE VYCPGSEKCP LVGSNVPWAF MQGEIATILA GDVKVKKERD P
Source	E.coli
Target Names	irf2bp2a
Protein Names	Recommended name: Interferon regulatory factor 2-binding protein 2-A Short name= IRF-2-binding protein 2-A Short name= IRF-2BP2-A
Expression Region	1-491
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.