



Recombinant Danio rerio Probable cytosolic iron-sulfur protein assembly protein ciao1 (ciao1)

Product Code	CSB-EP764705DIL-B
Abbreviation	ciao1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P0D9
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MKGALELQHR VSAHPDSRCW YVAWNPAGTT LATCGGDRAI RIWGKEGDSW ECKCVLSDGH QRTVRKVAWS PCGKYLASAS FDATTCIWKK TDEDFECLTV LEGHENEVKC VAWAPSGSLI ATCSRDKSVW IWEVDEEDEY ECLSVVNSHT QDVKHVVWHP TQELLASASY DNKICIYKEE DDDWECRATL EGHESTVWSL TFDPEGRRLA SCSDDRTVKI WKESTTGDGS SDESWKCICT LSGFHGRTIY DIAWCRLTGA LATACGDDGV RVFSEDPTAD PEQPIFALSA HVPKAHNQDV NCVSWNPKEA GLLATCSDNG EFAIWKYNSA
Source	E.coli
Target Names	ciao1
Protein Names	Recommended name: Probable cytosolic iron-sulfur protein assembly protein ciao1 Alternative name(s): WD repeat-containing protein 39
Expression Region	1-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.