



# Recombinant Rat Tyrosine-protein kinase Fgr (Fgr)

<b>Product Code</b>	CSB-EP764740RA
<b>Abbreviation</b>	Fgr
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6P6U0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	GCVFCKKLE PPKEDVGLG GDFRSQGAEE RYYPDPTQGR SSSISPQPIS PAFLNVGNIR SVSGTGVTIF VALYDYEART GDDLTFKGE KFHILNNTTEY DWWEARSLSS GRTGYVPSNY VAPVDSIQAE EWYFGKISRK DAERQLLSDG NPQGAFLIRE SETTKGAYSL SIRDWDQNRG DHIKHYKIRK LDMGGYYITT RAQFESVQDL VRHYMEVNDG LCYLLTAPCM VMKPQTLGLA KDAWEIDRNS IALDRRLGTG CFGDVWLGTW NCSTKVAVKT LKPGTMSPKA FLEEAIQIMKL LRHDKLVQLY AVVSEEPYI VTEFMCYGSL LDFLKDRKGH NLMLPNLVDM AAQVAEGMAY MERMNYIHRD LRAANILVGE HLIKIDAFG LARLIVDDEY NPQQGTFKPI KWTAPEAALF GRFTVKSDVW SFGILLTELI TKGRVPYPGM NNREVLEQVE HGYHMPCPPG CPVSLYEVME QTWRLDPEER PTFEYLQSFL EDYFTSTEPQ YQPGDQT
<b>Source</b>	E.coli
<b>Target Names</b>	Fgr
<b>Protein Names</b>	Recommended name: Tyrosine-protein kinase Fgr EC= 2.7.10.2 Alternative name(s): Proto-oncogene c-Fgr p55-Fgr
<b>Expression Region</b>	2-517
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.