



Recombinant Mouse Protein deltex-4 (Dtx4)

Product Code	CSB-BP764784MO
Abbreviation	Dtx4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6PDK8
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MLLASAVVW EWLNEHGRWR PYSPAVSHHI EAVVRAGPRA GGSVVLGQVD SRLAPYIIDL QSMNQFRQDT GTLRPVRRNY YDPSSAPGKG VWWEWENDNG SWTPYDMEVG ITIQYAYEKQ HPWIDLTSIG FSYIIDFSTM GQINRQTQRQ RRVRRRLDLI YPMVTGTMPK TQSWPVSPGP ATSSPAPPCS CPQCVLVMSV KAAVVHGGTG PPAVRKNMAL SGVGKLPQPP GPGAKPLDTT GTIRGPGKTA PSQVIRQVS NAPAGATVGS PASPQGSNRK TGRVALATLN RSNLQRLAIA QSRVLIASGV PTPVPKLNNG SSPVNPALAG ITGILMSAAG LPVCLTRPPK LVLHPPPVS SKSEIKSIPGVS NTSRKTTKKQ AKKGKTPEEV LKKYLQKVRH PPEEDCTICM ERLTAPSGYK GPQPTVKPDL VGKLSRCGHI YHIYCLVAMY NNGNKDGS LQ CPTCKTIYGV KTGTQPPGKM EYHLIPHS L P GHPDCKTIRI IYSIPPGIQG PEHPNPGKSF SARGFPRHCY LPDSEKGRKV LKLLLVAWDR RLIFAIGTSS TTGESDVIW NEVHHKTEFG SNLTGHGYPD ANYLDNVLAELAAQGIS EDS TSHEKD</p>
Source	Baculovirus
Target Names	Dtx4
Protein Names	Recommended name: Protein deltex-4 Short name= Deltex4
Expression Region	1-616
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.