



Recombinant Rat Protein BEX3 (Ngfrap1)

Product Code	CSB-MP764785RA
Abbreviation	Ngfrap1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6PDU5
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MANIQENEE MEQLQNGQE DRPVGGEGH QPAANNNHN HNHNHNHNHN HNHHRRGQAR RLAPNFRWAI PNRQMNDGLG GDGDDMEMFM EEMREIRKRL RELQLRNCLR ILMGELSNHH DHHDEFCLMP
Source	Mammalian cell
Target Names	Bex3
Protein Names	Recommended name: Protein BEX3 Alternative name(s): Brain-expressed X-linked protein 3 homolog Short name= rBex3 Nerve growth factor receptor-associated protein 1 p75NTR-associated cell death executor
Expression Region	1-130
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.