



Recombinant Human Serine protease FAM111B (FAM111B)

Product Code	CSB-BP764885HU
Abbreviation	FAM111B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6SJ93
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MNSMKTEENK SFSAMEDDQR TRPEVSKDTV MKQTHADTPV DHCLSGIRKC SSTFKLKSEV NKHETALEMQ NPNNLNNKECC FTFTLNGNSR KLDRSVFTAY GKPSESIYSA LSANDYFSER IKNQFNKNII VYEEKTIDGH INLGMPKCL PSDSHFKITF GQRKSSKEDG HILRQCENPN MECILFHVVA IGRTRKKIVK INELHEKGSK LCIYALKGET IEGALCKDGR FRSDIGEFEW KLKEGHKKIY GKQSMVDEVS GKVLEMDISK KKALQQKDIH KIKIQNESAT DEINHQSLIQ SKKKVHKPKK DGETKDVEHS REQILPPQDL SHYIKDKTRQ TIPRIRNYF CSLPRKYRQI NSQVRRRPHL GRRYAINLDV QKEAINLLKN YQTLNEAIMH QYPNFKEEAQ WVRKYFREEQ KRMNLSPAKQ FNIYKKDFGK MTANSVSVAT CEQLTYYSKS VGFMQWDNNG NTGNATCFVF NGGYIFTCRH VVHLMVGKNT HPSLWPDIIIS KCAKVTFTYT EFCPTPDNWF SIEPWLVKSN ENLDYAILKL KENGNAFPPG LWRQISPQPS TGLIYLIGHP EGQIKKIDGC TVIPLNERLK KYPNDCQDGL VDLYDTSNV YCMFTQRSFL SEVWNTHTLS YDTCFSDGSS GSPVFNASGK LVALHTFGLF YQRGFNVHAL IIEFGYSMDSI LCDIKKTNES LYKSLNDEKL ETYDEEKGKQ ESSLQDHQIE P MEC
Source	Baculovirus
Target Names	FAM111B
Protein Names	Recommended name: Protein FAM111B Alternative name(s): Cancer-associated nucleoprotein
Expression Region	1-734
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.