



Recombinant Danio rerio Tetratricopeptide repeat protein 35 (ttc35)

Product Code	CSB-MP764901DIL
Abbreviation	ttc35
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6TGY8
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSSISELYDV TWEEMRDCLR KWREENYRNS EQIVDVGEEL INEHASKLGD DIWIIYEQVM IAALDCSRDD LAWSCLQELK RQFPDSHRVK RLAGMRLEAL ERYDDANKLY DSILQDDPTN TAARKRKICI LRAQGKSSEA IRELNEYLEQ FVGDQEAWEH LSELYINEHD YAKAAFCLEE LMMTNPHNHL YCEQYAEVKY TQGGLENLEL ARKYFAQALK LNNRNMRALF GLYMSASHIA ASPKVSARKV KDNVKYAAWA AAQISRAYQM AGRGKKNKC SVKAVEEMLE SLQITMS
Source	Mammalian cell
Target Names	emc2
Protein Names	Recommended name: Tetratricopeptide repeat protein 35 Short name= TPR repeat protein 35
Expression Region	1-297
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.