



Recombinant Human Ran-binding protein 10 (RANBP10)

Product Code	CSB-EP764955HU-B
Abbreviation	RANBP10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6VN20
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AAATADPGA GNPQPGDSSG GGAGGGLPSP GEQELSRRLQ RLYPAVNQQE TPLPRSWSPK DKYNYIGLSQ GNLRVHYKGH GKNHKDAASV RATHPIPAAC GIYYFEVKIV SKGRDGYMGI GLSAQGVNMN RLPGWDKHSY GYHGDDGHSF CSSGTGQPYG PTFTTGDVIG CCVNLINGTC FYTKNGHSLG IAFTDLPANL YPTVGLQTPG EIVDANFGQQ PFLFDIEDYM REWRKAVQGT VHCFPISARL GEWQAVLQNM VSSYLVHHGY CATATAFARM TETPIQEEQA SIKNRQKIQK LVLEGRVGEA IETTQRFYPG LLEHNPMLL MLKCRQFVEM VNGTDSEVRS LSSRSPKSQD SYPGSPSLSP RHGPSSSHMH NTGADSPSCS NGVASTKSKQ NHSKYPAPSS SSSSSSSSSS SSPSSVNYSE SNSTDSTKSQ HHSSTSNQET SDSEMEMEAE HYPNGVLGSM STRIVNGAYK HEDLQTESS MDDRHPRRQL CGGNQAATER IILFGRELQA LSEQLGREYG KNLAHTEMLQ DAFSLAYS D PWSCPVGQQL DPIQREPVCA ALNSAILESQ NLPKQPPLML ALGQASECLR LMARAGLGSC SFARVDDYLH
Source	E.coli
Target Names	RANBP10
Protein Names	Recommended name: Ran-binding protein 10 Short name= RanBP10
Expression Region	2-620
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.