



# Recombinant Human Ran-binding protein 10 (RANBP10)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-YP764955HU   |
| <b>Abbreviation</b>      | RANBP10  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q6VN20   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Homo sapiens (Human)   |
| <b>Purity</b>            | ≥85% (SDS-PAGE)  |
| <b>Sequence</b>          | AAATADPGA GNPQPGDSSG GGAGGGLPSP GEQELSRRLQ RLYPAVNQQE<br>TPLPRSWSPK DKYNYIGLSQ GNLRVHYKGH GKNHKDAASV RATHPIPAAC<br>GIYYFEVKIV SKGRDGYMGI GLSAQGVNMN RLPGWDKHSY GYHGDDGHSF<br>CSSGTGQPYG PTFTTGDVIG CCVNLINGTC FYTKNGHSLG IAFTDLPANL<br>YPTVGLQTPG EIVDANFGQQ PFLFDIEDYM REWRKAVQGT VHCFPISARL<br>GEWQAVLQNM VSSYLVHHGY CATATAFARM TETPIQEEQA SIKNRQKIQK<br>LVLEGRVGEA IETTQRFYPG LLEHNPMLL MLKCRQFVEM VNGTDSEVRS<br>LSSRSPKSQD SYPGSPSLSP RHGPSSSHMH NTGADSPSCS NGVASTKSKQ<br>NHSKYPAPSS SSSSSSSSSS SSPSSVNYSE SNSTDSTKSQ HHSSTSNQET<br>SDSEMEMEAE HYPNGVLGSM STRIVNGAYK HEDLQTESS<br>MDDRHPRRQL CGGNQAATER IILFGRELQA LSEQLGREYG KNLAHTEMLQ<br>DAFSLAYS D PWSCPVGQQL DPIQREPVCA ALNSAILESQ NLPKQPPLML<br>ALGQASECLR LMARAGLGSC SFARVDDYLH |
| <b>Source</b>            | Yeast  |
| <b>Target Names</b>      | RANBP10  |
| <b>Protein Names</b>     | Recommended name: Ran-binding protein 10 Short name= RanBP10   |
| <b>Expression Region</b> | 2-620  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | Full Length of Mature Protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.