



Recombinant Danio rerio Myogenic factor 6 (myf6)

Product Code	CSB-EP764956DIL
Abbreviation	myf6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6VNZ9
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MMDLFETNAY FFNDLRYLEG DHGTLDMPGV SPLYEGNDSP LSPGQDPVPS ETGCESSGEE HVLAPPGLQA HCEGQCLMWA CKICKRKSAP TDRRKAATLR ERRRLKKINE AFDALKKKT V PNPQR LPKV EILRSAINYI EKLQDLLHSL DEQEQSNDTD PYTYNLKENH VTPSEYHWKK TCQSWQENPD HSSSQMAGHR EGAVLESSES SSLRRLSSIV DSISTEEPKA RCPSQISEK
Source	E.coli
Target Names	myf6
Protein Names	Recommended name: Myogenic factor 6 Short name= Myf-6 Alternative name(s): Muscle-specific regulatory factor 4
Expression Region	1-239
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.