



Recombinant Human Nasal embryonic luteinizing hormone-releasing hormone factor (NELF)

Product Code	CSB-EP764977HU-B
Abbreviation	NELF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6X4W1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	GAAASRRRA LRSEAMSSVA AKVRAARAFG EYLSQSHPEN RNGADHLLAD AYSGHDGSPE MQPAPQNKRR LSLVSNNGCYE GSLSEEPSIR KPAGEGPQPR VYTISGEPAL LPSPEAEAIE LAVVKGRRQR HPHHHSQPLR ASPGGSREDV SRPCQSWAGS RQGSKECPGC AQLAPGPTPR AFGLDQPPLP ETSGRRKKLE RMYSVDRVSD DIPIRTWFPK ENLFSFQTAT TTMQAISVFR GYAERKRRKR ENDSASVIQR NFRKHLRMVG SRRVKAQTFA ERRERSFSRS WSDPTPMKAD TSHDSRDSSD LQSSHCTLDE AFEDLDWDTE KGLEAVACDT EGFVPPKVML ISSKVPKAEY IPTIIRDDP SIIPILYDHE HATFEDILEE IERKLNVYHK GAKIWKMLIF CQGGPGHLYL LKNKVATFAK VEKEEDMIHF WKRLSRLMSK VNPEPNVIHI MGCYILGNPN GEKLFQNLRT LMPYRVTFE SPLELSAQGK QMIETYDFDR LYRLWKSQRH SKLLDFDDVL
Source	E.coli
Target Names	NSMF
Protein Names	Recommended name: Nasal embryonic luteinizing hormone-releasing hormone factor Short name= Nasal embryonic LHRH factor
Expression Region	2-530
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.