



# Recombinant Human Nasal embryonic luteinizing hormone-releasing hormone factor (NELF)

<b>Product Code</b>	CSB-EP764977HU
<b>Abbreviation</b>	NELF
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6X4W1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	GAAASRRRA LRSEAMSSVA AKVRAARAFG EYLSQSHPEN RNGADHLLAD AYSGHDGSPE MQPAPQNKRR LSLVSNHCYE GSLSEEPSIR KPAGEGPQPR VYTISGEPAL LPSPEAEAIE LAVVKGRRQR HPHHHSQPLR ASPGGSREDV SRPCQSWAGS RQGSKECPGC AQLAPGPTPR AFGLDQPPLP ETSGRRKKLE RMYSVDRVSD DIPIRTWFPK ENLFSFQTAT TTMQAISVFR GYAERKRRKR ENDSASVIQR NFRKHLRMVG SRRVKAQTFA ERRERSFSRS WSDPTPMKAD TSHDSRDSSD LQSSHCTLDE AFEDLDWDTE KGLEAVACDT EGFVPPKVML ISSKVPKAEY IPTIIRDDP SIIPILYDHE HATFEDILEE IERKLNVIYHK GAKIWKMLIF CQGGPGHLYL LKNKVATFAK VEKEEDMIHF WKRLSRLMSK VNPEPNVIHI MGCYILGNPN GEKLFQNLRT LMPYRVTFE SPLELSAQGK QMIETYDFDR LYRLWKSQRH SKLLDFDDVL
<b>Source</b>	E.coli
<b>Target Names</b>	NSMF
<b>Protein Names</b>	Recommended name: Nasal embryonic luteinizing hormone-releasing hormone factor Short name= Nasal embryonic LHRH factor
<b>Expression Region</b>	2-530
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.