



# Recombinant *Oryza sativa* subsp. japonica Putative B3 domain-containing protein Os08g0325100 (Os08g0325100, LOC\_Os08g23570)

<b>Product Code</b>	CSB-BP765037OFG
<b>Abbreviation</b>	Os08g0325100, LOC_Os08g23570
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6Z0D2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGTWCERCRR RDEQDYRNLD DCQKHFLLLM MGDFQHEIRG EQIVEQLTIP KEFVQRLKGD IPEEIQLETH NRNSYTVRVD KSQEKVIFAA GWAQFVKTFD LHMGDMMFR FKGNSQFDVI IFDQVGREKV CSVAVDDYLD PNVQEGRTDA TETLNSSRAH SQDDYLDPNV QEGRTNATET LNSSRAHSQP MPMQTPATET LNSSRAHSQD MPMQSPATET LNSSRAHPQP MPMQLPTETV NHFHAPHYPM QMPIENMALS RTQAMPTQM SPPTYRWTQV QRDNLRYSLP SEDQGCRVGV IPDPIGRRT KLNPVQEKVV NFKIQHIHSE IPIFVAVIKR SNVSGVLSTL SVAKRYVDEY LGGERFISLS RLGGKKGWIRL LLVGVGVAQG W
<b>Source</b>	Baculovirus
<b>Target Names</b>	Os08g0325134
<b>Protein Names</b>	Recommended name: Putative B3 domain-containing protein Os08g0325100
<b>Expression Region</b>	1-391
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.