



Recombinant *Oryza sativa* subsp. japonica Putative B3 domain-containing protein Os08g0325100 (Os08g0325100, LOC_Os08g23570)

Product Code	CSB-MP765037OFG
Abbreviation	Os08g0325100, LOC_Os08g23570
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6Z0D2
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MGTWCERCRR RDEQDYRNLD DCQKHFLLLM MGDFQHEIRG EQIVEQLTIP KEFVQRLKGD IPEEIQLETH NRNSYTVRVD KSQEKVIFAA GWAQFVKTFD LHMGDMMFR FKGNSQFDVI IFDQVGREKV CSVAVDDYLD PNVQEGRTDA TETLNSSRAH SQDDYLDPNV QEGRTNATET LNSSRAHSQP MPMQTPATET LNSSRAHSQD MPMQSPATET LNSSRAHPQP MPMQLPTETV NHFHAPHYPM QMPIENMALS RTQAMPTQM SPPTYRWTQV QRDNLRYSLP SEDQGCRVGV IPDPIGRRT KLNPVQEKVV NFKIQHIHSE IPIFVAVIKR SNVSGVLSTL SVAKRYVDEY LGGERFISLS RLGGKKGWIRL LLVGVGVAQG W
Source	Mammalian cell
Target Names	Os08g0325134
Protein Names	Recommended name: Putative B3 domain-containing protein Os08g0325100
Expression Region	1-391
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.