



Recombinant *Oryza sativa* subsp. *japonica* Rac-like GTP-binding protein 7 (RAC7)

Product Code	CSB-MP765046OFG
Abbreviation	RAC7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6Z7L8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MSTARFIKCV TVGDGAVGKT CMLISYTSNT FPTDYVPTVF DNFSANVVVD GNTVNLGLWD TAGQEDYNRL RPLSYRGADV FLLAFSLISK ASYENIHKKW IPELRHYAPN VPIVLVGTKL DLREDKQFFL DHPGLAPIST AQGEELKRMI GAAAYIECSS KTQQNVKSVF DSAIKVVLCPCPKPKKNTRK QRSC
Source	Mammalian cell
Target Names	RAC7
Protein Names	Recommended name: Rac-like GTP-binding protein 7 Alternative name(s): GTPase protein ROP5 OsRac7
Expression Region	1-194
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.